

**LAB ALERT - ANNOUNCEMENT**

***CURRENT NAME: HISTOPLASMA ANTIGEN DETECTION EIA, URINE***

***UPDATE NAME TO: HISTOPLASMA GALACTOMANNAN ANTIGEN QUANTITATIVE BY EIA, URINE***

**NOTIFICATION DATE: 1/2/2014**

**EFFECTIVE DATE: 1/6/2014**

<b>TEST CODE:</b>	UHSAG
<b>INTERFACE CODE:</b>	1001484
<b>METHODOLOGY:</b>	Quantitative Enzyme Immunoassay
<b>ANNOUNCEMENT:</b>	<p>The test has improved sensitivity and specificity for Histoplasma antigen detection. Previously it was reported in EIA units. The test is now quantitative and units are ng/mL. Previously detection cutoff was 2.0 EIA units.</p> <p>Updated reporting below:</p> <p>Reference Interval: Not Detected</p> <p>Interpretation Result: Detected or Not Detected</p> <p>If Antigen is Detected, result values will be reported as below:</p> <p>The quantitative range of this assay is 3.2 – 20.0 ng/mL and values will be reported within this range.</p> <p>Detection of Antigen outside of the linear range will not be reported as quantified values since results &lt;3.2 ng/mL or &gt; 20 ng/mL detection cannot be accurately quantified.</p>
<b>PERFORMING LAB:</b>	ARUP Laboratories
<b>ORDER/CONTACT INFORMATION:</b>	Thomas P. Lohmann, MD, Chief Medical Officer, 972-966-7135 or ClearPoint Client Services: 972-966-7700/Fax: 972-966-7799



Thomas P. Lohmann, MD  
 Chief Medical Officer