
LAB ALERT – INACTIVATION OF OLD TEST & REPLACEMENT WITH CURRENT TEST**NOTIFICATION DATE: 1/29/2014****EFFECTIVE DATE: 2/18/2014****REPLACEMENT TEST: BCR-ABL1 MUTATION ANALYSIS BY NEXT GENERATION SEQUENCING****OLD TEST: BCR-ABL1, T315I MUTATION DETECTION, QUANTITATIVE****REPLACEMENT TEST CODE: BCRNGS****REPLACEMENT INTERFACE CODE: 1004211****OLD TEST CODES TO BE DEACTIVATED: BCRT13****OLD INTERFACE CODES TO BE DEACTIVATED: 1004433****REASON FOR CHANGE: NEW MASSIVE PARALLEL SEQUENCING METHODOLOGY IS REPLACING THE POLYMERASE CHAIN REACTION METHODOLOGY****METHODOLOGY: MASSIVE PARALLEL SEQUENCING****PERFORMED: VARIES; REPORTED 10-12 DAYS****COLLECTION: Lavender (EDTA) Whole Blood or Bone Marrow (EDTA)****SPECIMEN PREPARATION: Transport 5 mL Whole Blood (Min: 1 mL) OR 3 mL Bone Marrow (Min: 1 mL)****STABILITY (FROM COLLECTION TO INITIATION OF TESTING): Ambient: 1 hour; Refrigerated: 48 hours; Frozen: Unacceptable****TRANSPORT: Refrigerated—Time Sensitive**

REFERENCE RANGE: Refer to Report

CPT CODE(S): 81403

PERFORMING LAB: ARUP Laboratories

ORDER/CONTACT INFORMATION: ClearPoint Client Services: 972-966-7700 | Fax: 972-966-7799
Thomas P. Lohmann, MD: 972-966-7135



Thomas P. Lohmann, MD
Chief Medical Officer